

Threshold testing in noise using a mobile, noise attenuating audiometric headset



**Odile Clavier, PhD
Jesse Norris, PhD
Brendan Flynn, PhD**

Create, LLC

**Sig Soli, PhD
AnnMarie Reebernacker**

House Clinic

Conflicts of Interest

- No conflicts of interest
- This research was supported by the National Institute On Deafness And Other Communication Disorders of the National Institutes of Health under Award Number R44DC012861. The content is solely the responsibility of the authors and does not necessarily represent the official views of the National Institutes of Health.


Increasing access to hearing healthcare

- Go mobile
- Put the audiometer in the earcups for calibrated stimuli
- Put the sound booth on the head
- Make it usable by untrained personnel




Integrated Audiometer Design

- 1. Deep earcup
- 2. Foam
- 3. Electronics
- 4. Speaker plate
- 5. Speaker and microphone
- 6. Grill
- 7. Fabric cover
- 8. Earseal and spacer



How much attenuation?



Adkisson MH. Ambient Noise Levels Present During School Hearing Screenings. University of Northern Colorado; 2012


Study Objectives

- Evaluate device performance in controlled noise environment
- Key question:
How much noise can device tolerate before we see threshold shifts?

Study Design


- 21 subjects (21 ears) with normal hearing (<20 dBHL)
- Study conducted at House Clinic in Los Angeles, CA
 - Subjects were seated in audiometric booths
 - Subjects completed 7 audiograms each
 - Twice in quiet (before and after noise conditions)
 - 47 dBA, 52 dBA, 57 dBA, 62 dBA, 67 dBA, brown noise
 - Noise condition order randomly selected for each subject

Noise Characteristics



Noise Characteristics


57 dBA condition



Threshold test


- Modified Hughson Westlake algorithm (automated)
- Step size of 2 dB
- Start level of 10 dB HL, maximum presentations = 35
- Traditionally: threshold determined as level where more than 50% of responses to ascending presentations are identical
- This can be too stringent with a 2dB step size – modified to: if three of five responses were within +/- 2 dB, the threshold is computed as the average of those three responses

Repeatability in quiet



Frequency (Hz)	250	500	1000	2000	3000	4000	6000	8000
Average difference (dB)	-0.9	-1.4	-0.8	0.4	0.2	-0.1	1.3	-0.9
Standard deviation (dB)	3.0	2.9	2.7	3.4	3.9	3.1	5.4	4.5

Results



OSHA criteria (industrial applications)

- Average thresholds at 2 – 4 kHz for each noise condition
- Compare each threshold to first measurement in quiet
- Shift = Differences of 10 dB or more
- ~~Average greater than 25 dB HL~~

Condition	% number of subjects
Quiet	0
47 dBA	0
52 dBA	0
57 dBA	4.8
62 dBA	28.6
67 dBA	71.4

Conclusions

- Normal thresholds can be reliably measured at ~55 dBA with only 5-10% of false positives according to several STS guidelines
- Screening at 20 dB HL could occur at even higher noise levels
- Adding high passive attenuation to the audiometric headset can greatly increase access to hearing testing across the world

